



**Training Title :** C-Arm

**Duration :** 2 Days

**Target group :** X-ray Technologists and Doctors

**Pre-requisites :** General Knowledge of X-ray, Anatomy

1st Day :

- 1- System Component
  - Acquisition Unit ;
  - Processing Unit ;
  - Options ;
  - Hands on.
- 2- System operation
  - System Movement ;
  - System ON/ OFF.
- 3- Software
  - Fluorospot Compact ;
  - Patient registration mode ;
  - Defining the examination settings ;
  - Acquiring images ;
  - Acquisition mode ;
  - Postprocessing subtask cards ;
  - Documentation ;
  - Data transfer.

2st Day :

- 1- Summary of the 1st day
- 2- Parameters adjustment
  - Image display ;
  - Examination properties ;
  - Image delete ;
  - DICOM properties ;
  - Site information.
- 3- Editing exams and organs programs
- 4- Practice with technicians and doctors